


# Naipo

## Instruction Manual Neck & Back MASSAGER

Neck Kneading / Back Swaying / Seat Vibration / Heating



Model: MGM-C11C

Shenzhen Pinlong E-Commerce Co., Ltd.  
Web: www.naipocare.com  
E-mail: support@naipocare.com  
Add: Huanan City Electronic Trading Building P05  
Room 116, Longgang District, Shenzhen, China



## Contents

English	01-03
Deutsch	04-06
Français	07-09
Español	10-12
Italiano	13-15

### TECHNICAL PARAMETERS

Name: NECK & BACK massager  
Model: MGM-C11C  
Voltage: DC12V  
Operating current: 4A  
Power: 48W  
Auto Timer: 15 minutes

### PRECAUTIONS

This unit is designed to be used by healthy adults. However, the following precautions should be taken when using the unit:

- Keep distance from humidity or fire.
- Before doing an exercise program, consult with a doctor to ensure that you are capable of conducting a regular routine without putting undue stress on your body.
- If you have never exercised before, had a surgery or suffer from any significant illness, including epilepsy, hemophilia, diabetes, cancer, hypertension, heart diseases, neuromuscular disorders or more, you should not conduct any exercise routines unless you have got written direction of your doctor.
- Do not use it while driving or sleeping.
- Do not use the unit by yourself.
- Do not use the unit under shower, bathtub, in case of any fire accident.
- Do not fit the unit with strap article.
- Do not use when you are outside.
- Use it only by children or pregnancy independently.
- Keep this unit out of reach of children. Do not connect it with any other power supply. Do not use this unit with any other electric equipment of the same timing, another electrical muscle stimulation device). Do not use this unit and any other machine at the same time. Please switch to "OFF" and isolate the power supply to stop using.
- Please use the specific path or adapter for the unit.
- Do not put the cable or plug on any hard article.
- This unit is not a medical or therapeutic device.
- The apparatus will not get hot on surface. Those sensitive to heat should pay special attention to it.

### FUNCTIONS

Kneading of neck, back, waist and thighs to release muscle tension and improve blood circulation. Double built-in vibrating motor of hip to relax the body and strengthen the metabolism. Unique up-down stroke is based on the theory of acupointology. The product was developed to be able to imitate human hands to knead and vibrate big area of your body with multi-massage heads in combination of kneading and vibration. Two pairs kneading heads with wadding function can knead, slide up and down to your body. There are three massage variations: full back, lower back and upper back massage, suitable to use in office, family and car.



Office Chair/Car Seat

-01-

-02-

-03-

-04-

-05-

-06-

-07-

-08-

-09-

-10-

-11-

-12-

-13-

-14-

-15-

-16-

-17-

-18-

-19-

-20-

-21-

-22-

-23-

-24-

-25-

-26-

-27-

-28-

-29-

-30-

-31-

-32-

-33-

-34-

-35-

-36-

-37-

-38-

-39-

-40-

-41-

-42-

-43-

-44-

-45-

-46-

-47-

-48-

-49-

-50-

-51-

-52-

-53-

-54-

-55-

-56-

-57-

-58-

-59-

-60-

-61-

-62-

-63-

-64-

-65-

-66-

-67-

-68-

-69-

-70-

-71-

-72-

-73-

-74-

-75-

-76-

-77-

-78-

-79-

-80-

-81-

-82-

-83-

-84-

-85-

-86-

-87-

-88-

-89-

-90-

-91-

-92-

-93-

-94-

-95-

-96-

-97-

-98-

-99-

-100-

-101-

-102-

-103-

-104-

-105-

-106-

-107-

-108-

-109-

-110-

-111-

-112-

-113-

-114-

-115-

-116-

-117-

-118-

-119-

-120-

-121-

-122-

-123-

-124-

-125-

-126-

-127-

-128-

-129-

-130-

-131-

-132-

-133-

-134-

-135-

-136-

-137-

-138-

-139-

-140-

-141-

-142-

-143-

-144-

-145-

-146-

-147-

-148-

-149-

-150-

-151-

-152-

-153-

-154-

-155-

-156-

-157-

-158-

-159-

-160-

-161-

-162-

-163-

-164-

-165-

-166-

-167-

-168-

-169-

-170-

-171-

-172-

-173-

-174-

-175-

-176-

-177-

-178-

-179-

-180-

-181-

-182-

-183-

-184-

-185-

-186-

-187-

-188-

-189-

-190-

-191-

-192-

-193-

-194-

-195-

-196-

-197-

-198-

-199-

-200-

-201-

-202-

-203-

-204-

-205-

-206-

-207-

-208-

-209-

-210-

-211-

-212-

-213-

-214-

-215-

-216-

-217-

-218-

-219-

-220-

-221-

-222-

-223-

-224-

-225-

-226-

-227-

-228-

-229-

-230-

-231-

-232-

-233-

-234-

-235-

-236-

-237-

-238-

-239-

-240-

-241-

-242-

-243-

-244-

-245-

-246-

-247-

-248-

-249-

-250-

-251-

-252-

-253-

-254-

-255-

-256-

-257-

-258-

-259-

-260-

-261-

-262-

-263-

-264-

-265-

-266-

-267-

-268-

-269-

-270-

-271-

-272-

-273-

-274-

-275-

-276-

-277-

-278-

-279-

-280-

-281-

-282-

-283-

-284-

-285-

-286-

-287-

-288-

-289-

-290-

-291-

-292-

-293-

-294-

-295-

-296-

-297-

-298-

-299-

-300-

-301-

-302-

-303-

-304-

-305-

-306-